

# **HGL Hydrophones**

The HGL Series hydrophones were designed to meet or exceed recommendations of section 3.3.2 of the AIUM Acoustic Output Measurement Standard (May 1998). They have an exceptionally flat sensitivity in a small and sturdy package. These hydrophones are excellent in-house standards for ultrasonic acoustic intensity measurements, and for general purpose field mapping.

#### **Features**

- High sensitivity
- Small effective aperture
- Broadband
- Solid construction
- Flawless integration with AH preamplifiers
- Flat (+/-3dB) 250 KHz to >> 20 MHz \*
- \* Use of the AG-20X0 preamplifier is required to maintain +/- 3 dB range from 20 to 40 MHz for the HGL-0085 and HGL-0200. For measurements above 40 MHz, the AH-20X0 preamplifier is advised.

#### **Technical Specifications**

	HGL-0085	HGL-0200	HGL-0400	HGL-1000
Frequency range (±3dB)	0.25 to 40 MHz		0.25 to 20 MHz	
Electrode aperture	85 μm	200 μm	400 μm	1000 μm
* EOC Nominal Sensitivity [nV/Pa]	8	45	160	510
Acceptance angle (-6dB at 5 MHz)	>150°	100°	30°	20°
Capacitance	30 pF			
Max. Operating Temperature	50 ℃			


<sup>\*</sup> EOC ("end of cable") is the open-circuit output sensitivity of the hydrophone. Calibration with an amplifier can be determined from the gain and input impedance of the amplifier.

Provided with traceable calibration 1-20 MHz at 50 KHz intervals. For other calibrations available visit our web site.




**HGL Hydrophone** 

### **Typical Sensitivity Plot**



### **Typical Directivity Plots**



## **Mechanical Specifications**

